

Functional potency assay predicts CAR-T effectiveness in tumor microenvironment


Yewei Xing*, Ningchun Liu*, Natalie Czerbya*, Evan Massi, Candy Garcia, Yelena Bronevetsky, Scott Wise*, James Lim

*Labcorp Early Development Laboratories Inc., Ann Arbor, Michigan; Xcellbio, San Francisco, California




Abstract: Predict CAR T potency in solid tumor microenvironment

BACKGROUND: Chimeric antigen receptor (CAR) T cell therapy holds great promise for the treatment of various cancers, including solid tumors. However, attempts to model the behavior and effectiveness of CAR-T cell therapies for blood cancers and solid tumors have been challenging due to the unique tumor microenvironments in which these cancer cells are found. The tumor microenvironment (TME) is often characterized by hypoxia, increased acidity, and high interstitial fluid pressures, allowing cancer cells to effectively evade immune surveillance. This immunosuppressive TME also contributes to CAR-T cell exhaustion, thereby limiting its antitumor activity and function.




METHODS: To address these concerns, we have developed a novel cell-based assay to measure CAR-T cell potency and cytotoxic function in three-dimensional (3D) *in vitro* cell culture system, human acute B cell lymphoblastic leukemia mouse model, and immunosuppressive tumor microenvironments. Utilizing the AVATAR system, we replicated the oxygen and interstitial fluid pressures found in the vasculature, the bone marrow and solid tumor microenvironments. Tumor cytotoxicity assays were conducted in these environments to measure cell exhaustion as analyzed by flow cytometry and electrical impedance. In addition, serial tumor challenge assays were performed to examine CAR-T potency and effectiveness in TME. RESULTS: Proof-of-concept experiments were performed using anti-receptor tyrosine kinase-like orphan receptor 1 (ROR1) CAR-T cells targeting the ovarian adenocarcinoma cell line, SKOV3. CD19 CAR-T were also used targeting the acute lymphoblastic leukemia cell line, NALM6. Defined ratios of effector T cells to tumor cells was assessed to model CAR-T potency *in vitro* and elevated CD19 CAR-T mediated cytotoxicity was confirmed with the increased ratio of effector T cells to tumor cells. Initial results from these screening experiments show significant decline in ROR1 CAR-T mediated cytotoxicity when performed under TME conditions. However, CD19 CAR-T showed effective cell killing under TME conditions. Interestingly, acclimating and expanding ROR1 CAR-T cells to high pressure and decreased oxygen culture (AVATAR) conditions improved potency levels and warrants further investigation. CONCLUSIONS: In summary, we observe CAR-T cells comprise the tumor cell killing ability in both *in vitro* and *in vivo* animal models. We also describe a physiologically relevant potency assay that incorporates hyperbaric and hypoxic incubation technology to predict the behavior of cell therapies in immunosuppressive tumor microenvironments.

Experimental design: Generating CAR-T cells under pressurized TME culture conditions to enhance potency in TME




Pressurized TME expansion condition enhances viral transduction, increasing the total number of CD19 CAR-T+ cells




Characteristics of CD19-CAR-T cells. Activated T cells were transfected with CD19-CAR lentivirus and cultured in standard or AVATAR incubators for 7 days before analyzed by flow cytometry. Cells grown in AVATAR system showed a 5% increase in transduction efficiency, together with a 20% increase in growth capacity.


TME expansion condition enhances CAR-T expansion rates, increasing the total number of CD19 CAR-T+ cells




Pressurized TME expansion condition enhances CD19-CAR-T survival during target tumor killing assay




Pressurized TME expansion condition maintains CD19 CAR-T potency during target tumor killing assay



Pressurized TME expansion condition enhances ROR1 CAR-T potency during serial tumor challenge assay



Pressurized TME expansion condition enhances ROR1 CAR-T potency under immunosuppressive TME screening conditions



Normoxia expanded CAR-T potency performance			AVATAR expanded CAR-T potency performance		
Target Killing Environment	Tumor Cytolysis % @ 24 hr	Tumor Cytolysis % @ 48 hr	Target Killing Environment	Tumor Cytolysis % @ 24 hr	Tumor Cytolysis % @ 48 hr
Normoxia	80.2	82.7	Normoxia	88.1 (+7.9%)	93.5 (+10.8%)
TME	69.1	57.7	TME	75.2 (+6.1%)	78.5 (+20.8%)

Scan the QR Code for a copy of this poster, and access to materials and methods section. Visit us at Xcellbio.com/avata-1 to find out how you can enhance your cell therapy potency in the solid tumor microenvironment.